Smart Use Smart Share RIOSensor





Your Clinic will be closer to SMART Digital Dentistry through the combination of RIOSensor and Windows tablet PC. It will improve efficiency of workflow and consultation.

Smart Use

Now, Share Intraoral X-ray images through any kind of smart phone, tablet, PC and laptop with RAYSCAN web providing necessary functions for most Intuitive user interface.



C€ **RIS500**

1113300		0120
Specifications	Size 1	Size 2
Outside dimensions Image area Pixel dimensions Sensor Technology Scintillator Number of pixels Image Resolution Connection Cable length	39x25mm 30x20mm 20x20um Advanced CMOS with optical fiber Csl. 1.5 megapixels 25 Ip/mm USB 2.0 high speed 2 meters	42x30mm 34x26mm 20x20um Advanced CMOS with optical fiber Csl. 2.2 megapixels 25 lp/mm USB 2.0 high speed 2 meters